



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz chromatography column with frit por. 0,
valve stopcock, sleeve NS 29/32, eff. length 600
mm, pouch 430 ml**

**€115.77
plus VAT &
Shipping**

Product Images



Description

Chromatography columns with fused-in frit (porosity:0), with valve stopcock (bore: 0-2.5 mm). Made of DURAN tube.

Additional Information

No.	LZ-1423
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3426
EAN	4058072311629
Transport temperature	Room temperature
diameter	30 mm
Material	Borosilicate glass
Volume	430 ml
Length	600 mm
Chromatography Columns Properties	with frit
Standard ground sleeve	29/32
Porosity	0
Type cock / valve	Valve tap

