



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz chromatography column with frit por. 0,
PTFE stopcock, core and sleeve NS 14/23, eff.
length 200 mm, 35 ml sticks**

**€122.77
plus VAT &
Shipping**

Product Images



Description

Chromatography columns with NS sleeve and NS core, with PTFE stopcock. With fused-in frit (porosity:0) With suction tube (thread GL 18), made of DURAN tube.

Additional Information

No.	LZ-1424
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3427
EAN	4058072311636
Transport temperature	Room temperature
diameter	15 mm
Material	Borosilicate glass
Volume	35 ml
Length	200 mm
Chromatography Columns Properties	with frit
Standard ground sleeve	14/23
Standard ground core	14/23
Porosity	0
Type cock / valve	PTFE - tap

