



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657




---

**Lenz chromatography column with frit por. 0,  
valve stopcock, core and sleeve NS 29/32, eff.  
length 400 mm, pou. 280 ml**

**€168.30  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

Chromatography columns with NS sleeve and NS core, with valve stopcock. With fused-in frit (porosity: 0). With suction tube (thread GL 18), made of DURAN tube.

## Additional Information

---

No.	LZ-1430
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3433
EAN	4058072311698
Transport temperature	Room temperature
diameter	30 mm
Material	Borosilicate glass
Volume	280 ml
Length	400 mm
Chromatography Columns Properties	with frit
Standard ground sleeve	29/32
Standard ground core	29/32
Porosity	0
Type cock / valve	Valve tap

