



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz chromatography column with frit por. 0,
valve stopcock, core and sleeve NS 29/32, eff.
length 600 mm, 430 ml capacity**

**€171.18
plus VAT &
Shipping**

Product Images



Description

Chromatography columns with NS sleeve and NS core, with valve stopcock. With fused-in frit (porosity: 0). With suction tube (thread GL 18), made of DURAN tube.

Additional Information

No.	LZ-1431
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3434
EAN	4058072311704
Transport temperature	Room temperature
diameter	30 mm
Material	Borosilicate glass
Volume	430 ml
Length	600 mm
Chromatography Columns Properties	with frit
Standard ground sleeve	29/32
Standard ground core	29/32
Porosity	0
Type cock / valve	Valve tap

