



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz chromatography column with frit Por. 0,
PTFE stopcock and sleeve NS 14/23, eff. length
200 mm, 35 ml,**

**€78.34
plus VAT &
Shipping**

Product Images



Description

High quality chromatography columns with standard ground joint, fused frit and PTFE stopcock. A cotton swab is not required. Available in various lengths and widths.

Additional Information

No.	LZ-1415
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3418
EAN	4058072311544
Transport temperature	Room temperature
diameter	15 mm
Material	Borosilicate glass
Volume	35 ml
Length	200 mm
Chromatography Columns Properties	with frit
Standard ground sleeve	14/23
Porosity	0
Type cock / valve	PTFE - tap

