



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657




**Lenz chromatography column with frit por. 0,  
PTFE stopcock and sleeve NS 29/32, eff. length  
400 mm, pou. 125 ml,**

**€93.56  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

High quality chromatography columns with standard ground joint, fused frit and PTFE stopcock. A cotton swab is not required. Available in various lengths and widths.

## Additional Information

---

No.	LZ-1416
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3419
EAN	4058072311551
Transport temperature	Room temperature
diameter	20 mm
Material	Borosilicate glass
Volume	125 ml
Length	400 mm
Chromatography Columns Properties	with frit
Standard ground sleeve	29/32
Porosity	0
Type cock / valve	PTFE - tap

