



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz chromatography column, 200x10 i. D., H 14/23, PTFE

€56.63
plus VAT &
Shipping

Product Images



Description

High-quality chromatography columns with punctures above the stopcock for placing a cotton plug, with NS sleeve, with PTFE plug and stopcock safety device. Made of borosilicate glass 3.3. Available in various lengths with different inner diameters.

Additional Information

No.	LZ-1402
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3405
EAN	4058072311414
Transport temperature	Room temperature
diameter	10 mm
Material	Borosilicate glass
Volume	15 ml
Length	200 mm
Standard ground sleeve	14/23
Type cock / valve	PTFE - tap

