



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz chromatography column, 400x20 i. D., H 29/32, PTFE

€68.25
plus VAT &
Shipping

Product Images



Description

High-quality chromatography columns with punctures above the stopcock for placing a cotton plug, with NS sleeve, with PTFE plug and stopcock safety device. Made of borosilicate glass 3.3. Available in various lengths with different inner diameters.

Additional Information

No.	LZ-1404
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3407
EAN	4058072311438
Transport temperature	Room temperature
diameter	20 mm
Material	Borosilicate glass
Volume	135 ml
Length	400 mm
Standard ground sleeve	29/32
Type cock / valve	PTFE - tap

