



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Lenz chromatography column, 800x40 i. D., H 29/32, PTFE

**€105.73**  
**plus VAT &**  
**Shipping**

## Product Images

---



## Description

---

High-quality chromatography columns with punctures above the stopcock for placing a cotton plug, with NS sleeve, with PTFE plug and stopcock safety device. Made of borosilicate glass 3.3. Available in various lengths with different inner diameters.

## Additional Information

---

No.	LZ-1406
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3409
EAN	4058072311452
Transport temperature	Room temperature
diameter	40 mm
Material	Borosilicate glass
Volume	1000 ml
Length	800 mm
Standard ground sleeve	29/32
Type cock / valve	PTFE - tap

