



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Desiccator plastic 27 cm, with tap, for plate 23 cm

€293.00
plus VAT &
Shipping

Product Images



Description

Desiccator made of polycarbonate, crystal clear, temperature resistant up to max 41°C.

Good stability due to large base area.

With tap for venting (housing made of PP, closure and sealing ring made of PTFE).

Tightly closing lid with silicone O-ring.

Holds vacuum (50 mbar) for 24 hours.

Additional Information

No.	2-5311
Manufacturer (Brand)	neoLab
EAN	4058072004576
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	273 mm
Height	311 mm
Height inside	197.0 mm
Desiccators properties	with tap with seal exsiccator plate available as accessory
Color	transparent (clear)
TBST MAX	41 °C
Material	Polycarbonate (PC)
Material gasket / septum	Silicone
Material tap	Polytetrafluoroethylene (PTFE)
Spare part	2-5315
Material closure/ lid	Polycarbonate (PC)
Position tap	on the lid
Type connection	Plug adapter. conical
for DM plate	230 mm

