



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Desiccator Glass 200 mm, NS 24/29 and valve

**€186.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The neoLab desiccator with standard ground joint 24/29 and valve is made of borosilicate glass 3.3 and is suitable for use under vacuum. Typical applications are the drying of samples and the storage of moisture-sensitive substances. Standardized according to DIN EN ISO 13130. Suitable perforated porcelain plates must be ordered separately.

## Additional Information

---

No.	E-1002
Manufacturer (Brand)	neoLab
EAN	4058072071080
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	270 mm
Height	300 mm
Desiccators properties	with tap with connection exsiccator plate available as accessory
Color	transparent (clear)
Material	Borosilicate glass
Material tap	Borosilicate glass Polytetrafluoroethylene (PTFE)
Standard ground	24/29
Material closure/ lid	Borosilicate glass
Position tap	on the lid
Type connection	Standard ground
for DM plate	190 mm
standardized according to DIN standard (DIN EN ISO)	DIN EN ISO 13130

