



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Desiccator Glass 300 mm, NS 24/29 and valve

€319.60
plus VAT &
Shipping

Product Images



Description

The neoLab desiccator with standard ground joint 24/29 and valve is made of borosilicate glass 3.3 and is suitable for use under vacuum. Typical applications are the drying of samples and the storage of moisture-sensitive substances. Standardized according to DIN EN ISO 13130. Suitable perforated porcelain plates must be ordered separately.

Additional Information

No.	E-1004
Manufacturer (Brand)	neoLab
EAN	4058072071103
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	380 mm
Height	424 mm
Desiccators properties	with tap with connection exsiccator plate available as accessory
Color	transparent (clear)
Material	Borosilicate glass
Material tap	Borosilicate glass Polytetrafluoroethylene (PTFE)
Material closure/ lid	Borosilicate glass
Standard ground	24/29
Position tap	on the lid
Type connection	Standard ground
for DM plate	280 mm
standardized according to DIN standard (DIN EN ISO)	DIN EN ISO 13130

