




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz ball flange, cup KF 25, thread M 16 x 1

€139.84
plus VAT &
Shipping



Description

Stainless steel ball flanges, threaded ladle

With KF cup, with thread

Additional Information

| | |
|---|------------------------|
| No. | LZ-1046 |
| Manufacturer (Brand) | Lenz |
| old neoLab article no. | B-3046 |
| EAN | 4058072307868 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| Color | silver |
| Material | Stainless steel |
| Form | straight |
| Connection General / Tempering jacket / Cooler / Vacuum | External thread M 16x1 |
| Spherical cut ball | KF 25 |

