



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Lenz ball flange, cup KF 25, thread M 30 x 1.5

**€165.58**  
**plus VAT &**  
**Shipping**



## Description

---

**Stainless steel ball flanges**, threaded ladle  
With KF cup, with thread

## Additional Information

---

No.	LZ-1047
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3047
EAN	4058072307875
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Material	Stainless steel
Form	straight
Connection General / Tempering jacket / Cooler / Vacuum	External thread M 30x1.5
Spherical cut ball	KF 25

