




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz gas sampling tube, capacity 150 ml, with sampling nozzle GL 14 and septum

€96.99
plus VAT &
Shipping

Product Images



Description

Gas sampling tubes according to DIN 12473-1

With sampling nozzle; with threaded connection GL 14 and septum, with two one-way stopcocks NS 14 (with threaded stopcock safety device).

Connections: Capillary tube (outer diameter 8 mm, inner diameter 2.7 mm), olive (outer diameter: 8 mm).

With labeling field

Additional Information

No.	LZ-1588
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3591
EAN	4058072313265
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	150 ml
Form	straight
Connection General / Tempering jacket / Cooler / Vacuum	GL 14
Type connection	Glass spigot
Type glass jar	Gas sample tube
Type valve	One way plug valve with NS plug

