



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz gas sampling tube, capacity 500 ml, with sampling nozzle GL 14 and septum

€116.17
plus VAT &
Shipping

Product Images



Description

Gas sampling tubes according to DIN 12473-1

With sampling nozzle; with threaded connection GL 14 and septum, with two one-way stopcocks NS 14 (with threaded stopcock safety device).

Connections: Capillary tube (outer diameter 8 mm, inner diameter 2.7 mm), olive (outer diameter: 8 mm).

With labeling field

Additional Information

| | |
|---|---------------------------------|
| No. | LZ-1590 |
| Manufacturer (Brand) | Lenz |
| old neoLab article no. | B-3593 |
| EAN | 4058072313289 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| Color | transparent (clear) |
| Material | Borosilicate glass |
| Volume | 500 ml |
| Form | straight |
| Connection General / Tempering jacket / Cooler / Vacuum | GL 14 |
| Type connection | Glass spigot |
| Type glass jar | Gas sample tube |
| Type valve | One way plug valve with NS plug |

