



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz gas sample tube, capacity. 500 ml, without
sampling tube, connections capillary tube and
olive**

**€102.40
plus VAT &
Shipping**

Product Images



Description

Gas sampling tubes according to DIN 12473-1

Without sampling port; with threaded connection GL 14 and septum, with two one-way stopcocks NS 14 (with threaded stopcock safety device).

Connections: Capillary tube (outer diameter 8 mm, inner diameter 2.7 mm), olive (outer diameter: 8 mm).

With labeling field

Additional Information

No.	LZ-1586
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3589
EAN	4058072313241
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	500 ml
Form	straight
Type connection	Glass spigot
Type glass jar	Gas sample tube
Type valve	One way plug valve with NS plug

