




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz gas inlet pipe, core NS 24/29, straight

€20.27
plus VAT &
Shipping

Product Images



Description

Gas inlet tubes, straight

Boiling capillaries, installation length: 200 mm. With core grind and thin extension tube (for pulling off at the burner).

Additional Information

No.	LZ-1003
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3003
EAN	4058072307431
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Form	straight
Length	200 mm
Standard ground core	24/29
Type cock / valve	Without

