




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz threaded pipe GL 32, 28 x 2.0 x 140 mm

€55.30
plus VAT &
Shipping

Product Images



Description

Threaded tubes for attaching with DIN thread (GL)

Additional Information

No.	LZ-0175
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2174
EAN	4058072299170
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass DURAN®
Length	140 mm
Connection General / Tempering jacket / Cooler / Vacuum	GL 32
standardized according to DIN standard (DIN)	DIN 12216

