



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz measuring tube without stopcock, volume 10 : 0.1 ml, zero point in tip

€113.05
plus VAT &
Shipping

Product Images

Description

Measuring tubes made of glass as spare parts for water determination apparatus according to Dean Stark
Zero point in the tip.

Additional Information

No.	LZ-1555
Manufacturer (Brand)	Lenz
old neolab article no.	B-3558
EAN	4058072312930
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass DURAN®
Volume	500 ml
Form	Y-shaped
Number of side necks	1
Quantity pistons	1
Standard ground sleeve	29/32
Standard ground core	29/32
Alignment side necks	Oblique 15°
Type cock / valve	Tap 1-Way
standardized according to DIN standard (DIN)	DIN 12420

