



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz dry pipe, U-shape, with NS stopcocks,
nominal length 125 mm, sleeve NS 14/23**

**€68.56
plus VAT &
Shipping**

Product Images



Description

Drying tubes according to DIN 12616, with replaceable NS stopcock plugs.
Side hose connections (AD: 6mm).
Made of DURAN® pipe.

Additional Information

No.	LZ-1306
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3309
EAN	4058072310455
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Form	U-shaped
Length	125 mm
Standard ground sleeve	14/23

