



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz cold trap with condensate drain, 2 connections NS 29/32, coolant content 500 ml

**€622.41
plus VAT &
Shipping**

Product Images



Description

Cold trap (cold traps) with condensate drain

For direct filling with coolants, with valve tap (with safety device) for condensate drain.

Double-walled design for low coolant consumption.

Apparatus connection: Sleeve NS 29/32

Vacuum pump connection: Core NS 29/32

Additional Information

No.	LZ-1311
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3314
EAN	4058072310509
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	500 ml
Standard ground sleeve	29/32
Standard ground core	29/32
Type cock / valve	Valve with PTFE piston

