



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz cold trap without condensate drain, 2 connections NS 29/32, coolant content 200 ml

€367.08
plus VAT &
Shipping

Product Images



Description

Cold traps (cold traps) without condensate drain For direct filling with refrigerants, double-walled design for low refrigerant consumption. Apparatus connection: Sleeve NS 29/32 Vacuum pump connection: Core NS 29/32

Additional Information

No.	LZ-1308
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3311
EAN	4058072310479
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	200 ml
Standard ground sleeve	29/32
Standard ground core	29/32
Type cock / valve	Without

