



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Lenz intermediate piece, cranked, sleeve lateral NS 14/23

**€65.16**  
**plus VAT &**  
**Shipping**

## Product Images

---



## Description

---

**Intermediate pieces**, cranked,  
with core and lateral sleeve NS 29/32, with sleeve NS 14/23 (for thermometer or boiling capillary).  
made from DURAN® tube

## Additional Information

---

No.	LZ-1280
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3283
EAN	4058072310196
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass DURAN®
Form	Y-shaped
Number of side necks	1
Standard ground sleeve	29/32
Standard ground core	29/32
Alignment side necks	Parallel
Standard ground for side necks Sleeve	14/23

