



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Lenz extraction apparatus for solids, extractor  
70 ml, with Dimroth condenser NS 45/40 with  
thread GL 14, flask 250 ml**

**€228.87  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

### Extraction apparatus for solids

For extraction from solids (warm)

#### Components:

Extraction attachment with extraction sleeve for solids

Dimroth condenser with thread GL 14

Round bottom flask

## Additional Information

---

|                               |  |
|-------------------------------|--|
| No.                           | LZ-1370                                      |
| Manufacturer (Brand)          | Lenz   |
| old neoLab article no.        | B-3373                                       |
| EAN                           | 4058072311094                                |
| Transport temperature         | Room temperature                             |
| autoclavable at 121°C         | Yes  |
| Color                         | transparent (clear)                          |
| Material                      | Borosilicate glass                           |
| Complete apparatus properties | with standard ground with cooler with piston |
| Quantity pistons              | 1  |
| Standard ground sleeve        | 45/40  |
| Standard ground core          | 29/32  |
| Type cooler                   | Dimroth cooler                               |
| Spare part                    | LZ-1373                                      |

