



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz extraction apparatus for solids, extractor
100 ml, with Dimroth condenser NS 45/40 with
thread GL 14, flask 250 ml**

**€232.94
plus VAT &
Shipping**

Product Images



Description

Extraction apparatus for solids

For extraction from solids (warm)

Components:

Extraction attachment with extraction sleeve for solids

Dimroth condenser with thread GL 14

Round bottom flask

Additional Information

No.	LZ-1371
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3374
EAN	4058072311100
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Complete apparatus properties	with standard ground with cooler with piston
Quantity pistons	1
Spare part	LZ-1374
Standard ground sleeve	45/40
Standard ground core	29/32
Type cooler	Dimroth cooler

