



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Lenz extraction apparatus for solids, extractor  
250 ml, with Dimroth condenser NS 45/40 with  
thread GL 14, flask 500 ml**

**€323.53  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

### Extraction apparatus for solids

For extraction from solids (warm)

#### Components:

Extraction attachment with extraction sleeve for solids

Dimroth condenser with thread GL 14

Round bottom flask

## Additional Information

---

No.	LZ-1372
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3375
EAN	4058072311117
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Complete apertures properties	with standard ground with cooler with piston
Quantity pistons	1
Standard ground sleeve	60/46
Standard ground core	29/32
Type cooler	Dimroth cooler
Spare part	LZ-1375

