



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz extraction apparatus for specific light
extraction agents, extractor 250 ml, cooler
connection NS 45/40, flask 500 ml**

**€497.09
plus VAT &
Shipping**

Product Images



Description

Extraction apparatus for specifically light extractants

For extraction from liquids containing specific light extractants.

Components:

Extraction attachment with insert filter (frit, porosity 0); with PTFE valve tap

Dimroth condenser with thread GL 14

Round bottom flask

Additional Information

No.	LZ-1380
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3383
EAN	4058072311193
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	250 ml
Complete apparatus properties	with standard ground with cooler
Standard ground sleeve	45/40
Standard ground core	29/32
Type cooler	Dimroth cooler
Spare part	LZ-1383

