



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz extraction apparatus for specifically heavy
extraction agents, extractor 100 ml, cooler
connection NS 45/40, flask 250 ml**

**€386.68
plus VAT &
Shipping**

Product Images



Description

Extraction apparatus for specifically heavy extractants
For extraction from liquids with specific heavy extractants.

Components:

Extraction attachment with insert filter (frit, porosity:0); with PTFE valve tap

Dimroth condenser with thread GL 14

Round bottom flask

Additional Information

No.	LZ-1388
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3391
EAN	4058072311278
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	100 ml
Complete apparatus properties	with standard ground with cooler with piston
Quantity pistons	1
Standard ground sleeve	45/40
Type cooler	Dimroth cooler
Spare part	LZ-1391

