



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz hot extraction apparatus acc. to Knöfler-
Böhm, extractor 30 ml, with Dimroth condenser
NS 29/32 with thread GL 14, flask 100 ml**

**€286.19
plus VAT &
Shipping**

Product Images



Description

Knöfler-Böhm **hot extraction apparatuses**

The condensate flows through the extraction material and continuously back into the flask.
The extracting agent is therefore in continuous circulation.

Components: Hot extraction attachment

Dimroth condenser with GL 14 thread

Round bottom flask

Additional Information

No.	LZ-1332
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3335
EAN	4058072310714
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	30 ml
Complete apparatus properties	with standard ground with cooler with piston
Quantity pistons	1
Standard ground sleeve	29/32
Standard ground core	29/32
Type cooler	Dimroth cooler
Spare part	LZ-1360

