



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz hot extraction apparatus acc. to Knöfler-
Böhm, extractor 250 ml, with Dimroth condenser
NS 45/40 with thread GL 14, flask 500 ml**

**€424.98
plus VAT &
Shipping**

Product Images



Description

Knöfler-Böhm **hot extraction apparatuses**

The condensate flows through the extraction material and continuously back into the flask. The extracting agent is therefore in continuous circulation.

Components: Hot extraction attachment

Dimroth condenser with GL 14 thread

Round bottom flask

Additional Information

No.	LZ-1334
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3337
EAN	4058072310738
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	250 ml
Complete apparatus properties	with standard ground with cooler with piston
Quantity pistons	1
Standard ground sleeve	45/40
Standard ground core	45/40
Type cooler	Dimroth cooler
Spare part	LZ-1361

