



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz hot extraction apparatus acc. to Knöfler-
Böhm, extractor 500 ml, with Dimroth
condenser NS 60/46 with thread GL 14, flask
1,000 ml**

**€590.87
plus VAT &
Shipping**

Product Images



Description

Knöfler-Böhm **hot extraction apparatuses**

The condensate flows through the extraction material and continuously back into the flask.
The extracting agent is therefore in continuous circulation.

Components: Hot extraction attachment

Dimroth condenser with GL 14 thread

Round bottom flask

Additional Information

No.	LZ-1335
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3338
EAN	4058072310745
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	500 ml
Complete apparatus properties	with standard ground with cooler with piston
Quantity pistons	1
Spare part	LZ-1362
Standard ground sleeve	60/46
Standard ground core	60/46
Type cooler	Dimroth cooler

