



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz extraction apparatus acc. to Soxhlet,
extractor 150 ml, with Allihn condenser NS
45/40, flask 250 ml**

**€267.73
plus VAT &
Shipping**

Product Images



Description

Soxhlet extraction apparatus with Allihn condenser

For extraction from solids with a suitable solvent.

This method allows effective separation with relatively small amounts of extractant.

Components: Extractor with piston connection NS 29/32

Allihn condenser with thread GL 14

Round bottom flask

Additional Information

No.	LZ-1324
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3327
EAN	4058072310639
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	150 ml
Complete apparatus properties	with standard ground with cooler with piston
Quantity pistons	1
Standard ground sleeve	45/40
Standard ground core	45/40
Type cooler	Allihn cooler
Spare part	LZ-1343 LZ-1351

