



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz extraction apparatus acc. to Soxhlet,
extractor 300 ml, with Dimroth condenser NS
60/46, flask 500 ml**

**€391.80
plus VAT &
Shipping**

Product Images



Description

Soxhlet extraction apparatus with Dimroth condenser

For extraction from solids with a suitable solvent.

This method allows effective separation with relatively small amounts of extractant.

Components: Extractor with piston connection NS 29/32

Dimroth condenser with thread GL 14

Round bottom flask

Additional Information

No.	LZ-1318
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3321
EAN	4058072310578
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	300 ml
Complete apparatus properties	with standard ground with cooler with piston
Standard ground sleeve	60/46
Standard ground core	60/46
Type cooler	Dimroth cooler
Spare part	LZ-1346

