



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz extraction apparatus acc. to Soxhlet with
flange, capacity 2,000 ml, flat ground joint lid LF
100, clamping closure LF 100, with Dimroth
condenser NS 45/40, flask 2,000 ml NS 29/32**

**€890.59
plus VAT &
Shipping**

Product Images



Description

Extraction apparatus according to Soxhlet, 2000 ml, with flange

The apparatus consists of:

Extraction attachment with flange LF 100

Plane ground cover LF 100; middle neck NS 45/40

Clamping cap LF 100

Dimroth condenser with core NS 45/40; threaded connections GL 14

Round bottom flask with NS 29/32; capacity: 2000 ml

Additional Information

No.	LZ-1327
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3330
EAN	4058072310660
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	2000 ml
Complete apparatus properties	with standard ground with cooler with piston
Quantity pistons	1
Spare part	LZ-1328
Standard ground sleeve	71/51
Standard ground core	71/51
Type cooler	Dimroth cooler

