



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz extraction apparatus according to
Thielepape, extractor 250 ml, with Dimroth
condenser NS 45/40, extraction attachment
with glass frit sleeve, flask 500 ml**

**€447.78
plus VAT &
Shipping**

Product Images



Description

Extraction apparatus according to Thielepape

For the extraction of solid substances

The apparatus consists of:

Extraction attachment with glass frit sleeve

Dimroth condenser with thread GL 14

Round bottom flask

Additional Information

No.	LZ-1331
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3334
EAN	4058072310707
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	250 ml
Complete apparatus properties	with standard ground with cooler with piston
Quantity pistons	1
Standard ground sleeve	45/40
Standard ground core	45/40
Type cooler	Dimroth cooler
Spare part	LZ-1356

