



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz hot extraction apparatus according to
Twisselmann, extractor 150 ml, cooler
connection NS 45/40, extraction attachment
with vacuum jacket, flask 250 ml**

**€668.47
plus VAT &
Shipping**

Product Images



Description

Hot extraction apparatuses according to Twisselmann

Components:

Hot extraction attachment with vacuum jacket

Twisselmann cooler (improved) with thread GL 14

Round bottom flask

Additional Information

No.	LZ-1339
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3342
EAN	4058072310783
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	150 ml
Complete apparatus properties	with standard ground with cooler with piston
Quantity pistons	1
Standard ground sleeve	45/40
Standard ground core	45/40
Type cooler	Twisselmann cooler
Spare part	LZ-1369 LZ-1366

