




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz reaction apparatus with descending distillation, reaction vessel with laboratory flange LF 150, capacity 10,000 ml, with temperature control jacket, complete glassware set, stirrer and stand

**€16,516.60
plus VAT &
Shipping**



Description

Reaction apparatuses with descending distillation

For discontinuous reaction, for stirring under normal pressure or vacuum, with the possibility of adding liquids and boiling under reflux. The descending cooler with intermediate and final receiver (suitable for vacuum operation) allows the separation of the solvent from the reaction product via one-way distillation after the end of the reaction.

High quality system components

Reaction vessels (6,10,20 or 30l)

with temperature control jacket (heatable via circulation thermostat), with bottom drain valve (PTFE, no dead volume)

Plane ground cover (LF 150 or LF 200)

1 middle neck NS 45/40, 2 side necks 45° NS 29/32, 1 side neck straight NS 14/23

Dosing funnel

Capacity: 1000 or 4000 ml with graduation, NS 29/32

Spiral cooler

NS 29/32, S 29, jacket length: 300 mm

Stirring closure and stirrer

with 16 mm stirrer shaft, stirrer: steel, PTFE-coated

Additional Information

No.	LZ-1934
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3936
EAN	4058072316716
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass steel
Volume	10 l
sterile	No
Wide	500 mm
Height	2000 mm
Length	1000 mm
Depth	1000 mm
Type closure/ lid	Flat ground cover LF 150
Complete apparatus properties	in set
Drive	Agitator with motor 230V and stepless regulation
Number of parts in set	48
Spare part	B-3210 LZ-1814
Stirrer	Stirrer with movable blades and coupling pins
Stirring caps	Stirring caps with double lip seal
Type cooler	Dimroth cooler
Type reaction vessel	Reaction vessel with temperature control jacket. Bottom drain valve
Type tripod	Pilcoated steel pipe. FL 150. complete with clamps and sockets
Type funnel	Dropping funnel. cylindrical. PTFE plug. Pressure equalization tube

