



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz water determination apparatus according to Dean Stark, 0-point above the tap, measuring tube with tap

€202.71
plus VAT &
Shipping

Product Images



Description

Water determination apparatus according to Dean Stark according to DIN 12420 (xylene method), volume 10:0.1 ml, with center-neck round-bottom flask 500 ml, ground joints: NS 29/32, measuring tube with or without stopcock, with Liebig cooler, jacket length 400 mm

zero point above stopcock

Additional Information

No.	LZ-1551
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3554
EAN	4058072312893
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass DURAN®
Volume	500 ml
Form	Y-shaped
Height	400 mm
Position zero point	above the cock
Standard ground sleeve	29/32
Standard ground core	29/32
Spare part	LZ-1554
Type connection	Standard grind hose olive
Type ground	Standard ground
Type valve	One way plug valve with NS plug

