



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz cooler according to Davies with olives,
sleeve/core NS 24/29, shell length 400 mm,
cooling surface 740 cm²**

**€139.07
plus VAT &
Shipping**

Product Images



Description

Standard ground joint reflux condenser for condensation of gases and vapors in chemical synthesis. These high-quality reflux condensers made of borosilicate glass 3.3 are nowadays indispensable in the chemical synthesis laboratory. Cooler according to Davies: With glass tubing olives. Suitable for tubing with an inner diameter of 9–12 mm.

Additional Information

No.	LZ-1127
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3129
EAN	4058072308674
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Length	0 400
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection cooling medium	Glass olive
Cooling surface	740 cm ²
Material olives	Borosilicate glass
Standard ground sleeve	24/29
Standard ground core	24/29
Type cooler	Davies cooler

