



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Lenz distillation bridge with Liebig condenser, vacuum advance and threaded connections GL 14 (cooling water) and GL 18 (vacuum), cores 14/23, length 160 mm, cooling surface 60 cm<sup>2</sup>**

**€119.15  
plus VAT &  
Shipping**

## Product Images



## Description

**Distilling bridges with Liebig cooler.** Drain pipe with NS core, threaded connections: GL 14 (cooling water), GL 18 (vacuum).

With ground joint sleeve NS 14/23 for distillation thermometer with installation length 55 mm.

Made of DURAN® pipe.

## Additional Information

No.	LZ-0971
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2971
EAN	4058072307110
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Form	T-shaped angled
Length	160 mm
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection General / Tempering jacket / Cooler / Vacuum	GL 14 GL 18
Standard ground sleeve	14/23
Standard ground core	14/23
Type cooler	Liebig cooler

