



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz distillation bridge with Liebig condenser, vacuum advance and threaded connections GL 14 (cooling water) and GL 18 (vacuum), cores 14/23, length 250 mm, cooling surface 100 cm²

**€130.65
plus VAT &
Shipping**

Product Images



Description

Distilling bridges with Liebig cooler. Drain pipe with NS core, threaded connections: GL 14 (cooling water), GL 18 (vacuum).

With ground joint sleeve NS 14/23 for distillation thermometer with installation length 55 mm.

Made of DURAN® pipe.

Additional Information

No.	LZ-0974
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2974
EAN	4058072307141
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Form	T-shaped angled
Length	250 mm
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection General / Tempering jacket / Cooler / Vacuum	GL 14 GL 18
Standard ground sleeve	14/23
Standard ground core	14/23
Type cooler	Liebig cooler

