



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz distillation bridge for short path
distillation, 2 cores NS 14/23, 1 sleeve NS 14/23
and 3 threaded connections GL 14**

**€110.02
plus VAT &
Shipping**

Product Images



Description

Short condensation path due to compact design of cooler and vacuum connection, with two NS cores 14/23, one NS sleeve 14/23 and three threaded connections GL 14.
Made of DURAN® tube.

Additional Information

No.	LZ-0994
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2994
EAN	4058072307349
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection General / Tempering jacket / Cooler / Vacuum	GL 14
Standard ground sleeve	14/23
Standard ground core	14/23
Type cooler	Liebig cooler

