



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz distillation bridge, cores NS 14/23, length 150 mm

€30.15
plus VAT &
Shipping

Product Images



Description

Without cooler, with two vertical NS cores.
Made of DURAN® tube.

Additional Information

No.	LZ-0950
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2950
EAN	4058072306908
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Form	angled
Length	150 mm
Cooler features	with standard ground joint
Material	Borosilicate glass
Standard ground core	14/23

