



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Lenz distillation bridge with Liebig condenser,  
cores NS 14/23, cooling water connection GL 14,  
length 250 mm, cooling surface 100 cm<sup>2</sup>**

**€104.72  
plus VAT &  
Shipping**

## Product Images



## Description

---

**Distilling bridges with Liebig cooler;** drain pipe with NS core, cooling water connection with DIN thread GL 14, with ground joint sleeve NS 14/23 for distillation thermometers with an installation length of 55 mm.  
Made of DURAN® pipe.

## Additional Information

---

No.	LZ-0966
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2966
EAN	4058072307066
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Form	T-shaped
Length	250 mm
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection General / Tempering jacket / Cooler / Vacuum	GL 14
Standard ground core	14/23
Type cooler	Liebig cooler

