



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz distillation bridge with ground socket NS 14/23, cores NS 19/26, length 200 mm

€53.73
plus VAT &
Shipping

Product Images



Description

Distilling bridge with two vertical NS cores; without cooler.

Ground joint sleeve 14/23 for distillation thermometers with an installation length of 55 mm.

Made of DURAN® tube.

Additional Information

No.	LZ-0959
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2959
EAN	4058072306991
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Form	T-shaped angled
Length	200 mm
Cooler features	with standard ground joint
Material	Borosilicate glass
Standard ground sleeve	14/23
Standard ground core	19/26

