



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz distillation bridge according to Claisen,
with GL connection, with Liebig cooler, cores
14/23, with cooling water connection GL 14,
length 160 mm, cooling surface 60 cm²**

**€107.07
plus VAT &
Shipping**

Product Images



Description

Distilling bridges with Liebig cooler (straight outlet) and attachment according to Claisen, cooling water connection with DIN thread GL 14, with two ground joint sleeves NS 14/23 for distillation thermometers with an installation length of 55 mm.

Made of DURAN® tube.

Additional Information

No.	LZ-0979
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2979
EAN	4058072307196
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Length	160 mm
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection General / Tempering jacket / Cooler / Vacuum	GL 14
Standard ground sleeve	14/23
Standard ground core	14/23
Type cooler	Liebig cooler

