



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz distillation bridge according to Claisen,
with Liebig condenser, cores NS 14/23, cooling
water connection olive, length 250 mm, cooling
surface 100 cm²**

**€98.57
plus VAT &
Shipping**

Product Images



Description

Distilling bridges with Liebig cooler (straight outlet) and attachment according to Claisen, cooling water connection with olive, with two ground joint sleeves NS 14/23 for distillation thermometers with an installation length of 55 mm. Made of DURAN® tube.

Additional Information

No.	LZ-0980
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2980
EAN	4058072307202
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Length	250 mm
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection General / Tempering jacket / Cooler / Vacuum	Olive
Standard ground sleeve	14/23
Standard ground core	14/23
Type cooler	Liebig cooler

