



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz distillation bridge acc. to Claisen, with vacuum advance and GL connection, with Liebig condenser, cores NS 29/32, cooling water connection GL 14, vacuum connection GL 18, length 250 mm, cooling surface 130 cm²

**€148.14
plus VAT &
Shipping**

Product Images



Description

Distilling bridges according to DIN 12594, with Liebig cooler and attachment according to Claisen, threaded connections: GL 14 (cooling water) and GL 18 (vacuum), with two ground joint sleeves NS 14/23 for distillation thermometers with an installation length of 55 mm.

Made of DURAN® tube.

Additional Information

No.	LZ-0991
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2991
EAN	4058072307318
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Length	250 mm
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection General / Tempering jacket / Cooler / Vacuum	GL 14 GL 18
Standard ground sleeve	14/23
Standard ground core	29/32
Type cooler	Liebig cooler

