



neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Lenz distillation bridge acc. to Claisen, with vacuum advance and GL connection, with Liebig condenser, cores NS 29/32, cooling water connection GL 14, vacuum connection GL 18, length 400 mm, cooling surface 210 cm<sup>2</sup>**

**€146.66  
plus VAT &  
Shipping**

## Product Images



## Description

---

**Distilling bridges** according to DIN 12594, with Liebig cooler and attachment according to Claisen, threaded connections: GL 14 (cooling water) and GL 18 (vacuum), with two ground joint sleeves NS 14/23 for distillation thermometers with an installation length of 55 mm.

Made of DURAN® tube.

## Additional Information

---

No.	LZ-0993
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2993
EAN	4058072307332
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Length	400 mm
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection General / Tempering jacket / Cooler / Vacuum	GL 14 GL 18
Standard ground sleeve	14/23
Standard ground core	29/32
Type cooler	Liebig cooler

